

Product Features

MANYFILL D VIAL

Filling/injection multipatient* system with connection for contrast media and saline solution in bottles and clamps

For all dual syringe CT scanner injectors

*CAUTION : The Manyfill range enables filling and injection in close-system for several patients, only if there is a patient line secufill protection downstream. Secufill must be changed for each patient (consult the instructions for use).



Technical data :

Reference	MANYFILL D VIAL
Reservoir	Yes
Luer	yes conform ISO 594
Length	150 cm + 20 cm
Filling/Contrast media: Yellow tube	Ø inner 3,0 mm Ø outer 4,1 mm
Filling/Saline solution : Blue tube	Ø inner 2,3 mm Ø outer 4,1 mm
Injection/Patient: Colorless tube	Ø inner 2,3 mm Ø outer 4,1 mm
Opening threshold	30 kPa (4,2 PSI)
Maximum pressure	2100 kPa (300 PSI)
Compatibility	all dual syringe CT injector
Material	PVC, PE, PP, ABS, PA, silicone
Council directive 93 / 42 / EEC	II a
Sterilization	Ethylene Oxide (EO)
Phthalate free	



The Manyfill D Vial has 2 one way valves.

Packaging : 15 Units

Pack content :

	Contrast media		Saline solution	
	Vial	Luer lock bag	Vial	Luer lock bag
MANYFILL D VIAL	1)		1)	

Operating conditions:



→ Injection way
→ Filling way

